

# himac APPLICATION

No. 60 April 1997


**Subject** Separation of plasmid DNA by neo-angle rotor

**Model** Preparative Micro-Ultracentrifuge CS-GX series  
Neo-angle rotor S120NT


Rapid separation of plasmid DNA using 2-step gradient centrifugation by neo-angle rotor S120NT (holding 2ml tube).

We've shown in our himac APPLICATION it takes 2h to separate plasmid DNA by neo-angle rotor S120NT <sup>1)</sup>. We tried this separation even shorter using 2-step gradient centrifugation. As a result we enabled this separation for 1 hour.

## Results



Before centrifugation



After centrifugation

## Conditions

### (1) Condition of ultracentrifugation

Rotor	Speed (rpm)	Time (h)	Temp. (°C)	Acc. mode	Dec. mode
S120NT neo-angle rotor	120,000	1.0	20	"5"	"7"

### (2) Used tube : 2PA Cone-Top tube\*

\*Registered trade mark of Seton Scientific Company.

### (3) Sample preparation

We used E.coli., JM109 holding plasmid pUC 19 DNA. The E.coli., JM109 was cultured over night and then crude DNA was isolated by alkalin-SDS method and other processes and dissolved TE buffer. We used this solution as a sample.

For each 2PA Cone-Top tube\*, we provided :

(Low density solution : Upper layer)

42% (w/w) CsCl/TE buffer (pH8.0) ; about 0.9ml

(High density solution : Lower layer)


Sample: 0.62ml

CsCl : 0.98g

Ethydium Bromide (10mg/ml) : 40 $\mu$ l

Polyethylen (10) octylphenylether (Triton X - 100) : 2 $\mu$ l

Add 0.9ml of low density solution into a 2PA Cone-Top tube\* using a syringe. And add about 0.9ml of high density solution gently from the tube bottom. (Use a syringe equipped to with a long needle : see below figure.) If it is not enough to fill the tube the tube should be filled with low density solution.



#### Reference

- 1) himac APPLICATION No. 59 (1997).

# HITACHI

#### (Export office)

Nissei Sangyo Co., Ltd.

Head Office :

1 - 24 - 14, Nishishinbashi, Minato-ku, Tokyo, 105 Japan

Tel : 81 - 3 - 3504 - 7281

Fax : 81 - 3 - 3504 - 7302

#### (Manufacturer)

Hitachi Koki Co., Ltd.

Scientific Instruments Div.

1060, Takeda, Hitachinaka - city, Ibaraki - pref., 312 Japan

Tel : 81 - 29 - 276 - 7384

Fax : 81 - 29 - 276 - 7475