

# himac APPLICATION

No. 62 April 1997

**Subject** Separation of chloroplast DNA

**Model** Preparative Micro-Ultracentrifuge CS-GX series  
Fixed-angle rotor S100AT6

Separation of chloroplast DNA from spinach by fixed-angle rotor S100AT6

It is known chloroplast and mitochondria have own DNA distinct from chromosomal DNA. Here we tried to separate spinach chloroplast DNA by ultracentrifugation. As a result it took over night (16 hours) to separate by CsCl-density gradient ultracentrifugation using preparative micro-ultracentrifuge CS-GX series and fixed-angle rotor S100AT6 (holding 5ml tube).

**Results**




Fig.1 Before centrifugation




Fig.2 After centrifugation

**Conditions**

(1) Condition of ultracentrifugation

| Rotor                     | Speed (rpm) | Time (h) | Temp. (°C) | Acc. mode | Dec. mode |
|---------------------------|-------------|----------|------------|-----------|-----------|
| S100AT6 fixed-angle rotor | 70,000      | 16       | 20         | "9"       | "7"       |

(2) Used tube : 5PA Re-seal tube\*

\*Registered trade mark of Seton Scientific Company.

### (3) Sample preparation<sup>1)</sup>

5ml of chloroplast suspension (separated as described <sup>2)</sup>)

↓  
← Add 5 volumes of buffer A.

↓  
centrifugation

↻ CF8DL refrigerated centrifuge, 2,500rpm, 15min, 4°C, 50ml volume tube.

↓  
pellet

↓  
← Wash the pellet for 2 or 3 times with 20-30ml of buffer A.

↓  
washed pellet

↓  
← Resuspend 2ml of buffer A and add pronase.

↓  
incubation

↓  
← Incubate 37°C, 2h and then room temp. 2min.

↓  
mix

↓  
← Mix gently for 15min.

↓  
incubation

↓  
← Incubate 4°C for 2-3 hours.

↓  
centrifugation

↻ CF8DL refrigerated centrifuge, 2,500rpm, 15min, 4°C.

↓  
supernatant (crude DNA solution)

↓  
ultracentrifugation

↻ Fixed-angle rotor S100AT6, 70,000rpm, 16h, 20°C.

↓  
chloroplast DNA

Buffer A : 0.35M Sorbitol, 50mM Tris-HCl (pH8.0), 25mM EDTA-Na<sub>2</sub>.

Buffer B : 5% (w/w) Sodium N-lauroyl Sarcosinate, 50mM Tris-HCl (pH8.0), 25mM EDTA Na<sub>2</sub>.

#### Reference

1) METHODS IN ENZYMOLOGY, 118, 167 - 186 (1986).

2) himac APPLICATION No.61 (1997).

# HITACHI

#### (Export office)

Nissei Sangyo Co., Ltd.

Head Office :

1 - 24 - 14, Nishishinbashi, Minato-ku, Tokyo, 105 Japan

Tel : 81 - 3 - 3504 - 7281

Fax : 81 - 3 - 3504 - 7302

#### (Manufacturer)

Hitachi Koki Co., Ltd.

Scientific Instruments Div.

1060, Takeda, Hitachinaka - city, Ibaraki - pref., 312 Japan

Tel : 81 - 29 - 276 - 7384

Fax : 81 - 29 - 276 - 7475