

himac APPLICATION

No. 76 June 1998

Subject Separation of total RNAs from homogenated cell line with fixed-angle rotor

Model S150AT fixed-angle rotor for preparative micro-ultracentrifuge CS150GX

Separation example of total RNAs from homogenated cell line using fixed-angle rotor only in 1.75 hours

Most of literature describe that a separation of total RNAs from homogenated cell line has been generally processed by overnight ultracentrifugation using swinging bucket rotor.¹⁾ On the contrary, Hitachi has been stated that separation time for total RNAs can be reduced up to 2.5 hours by using micro-ultracentrifuge and its RP120AT fixed-angle rotor.²⁾ Moreover, we have studied to reduce separation time much more by using newly developed micro-ultracentrifuge, model CS150GX of which maximum speed is 150,000rpm. As a result of the study, we could know that a speed of 150,000rpm can separate total RNAs only in 1.75 hours. By the experimental result of electrophoresis and RT-PCR, obtained total RNAs in such a rapid separation can be equivalent to that obtained by the former centrifugation.

1.Results of experiment




Fig.1 After centrifugation
(Finally, 1 μ g/ml of EtBr added)

Fig.2 Result of electrophoresis

1:Separation with S100AT6
2:Separation with S150AT




Fig.3 Electrophoresis of DNA amplified by RT-PCR

2.Contents of experiment

(1) Conditions of centrifugation

Rotor	Speed	Running time	Temperature	Acc.	Dec.
S150AT fixed-angle	150,000rpm	1.75 hours	20°C	5	7

(2) Used tube: 2PA Cone-Top tube* (* Registered trade mark of Seton Scientific Company.)

(3) Methods of experiment

Cell pellet(5-10 x 10⁶ cells/tube)



4M GTC 1-1.5ml, Vortex



Stack on 5.3M CsCl sol* 0.8ml



Ultracentrifugation



Dissolve into buffer and precipitate by ethanol

*Separation of total RNAs may not be completed due to crystallization during centrifugation when 5.7M CsCl is used.

(References)

- 1) Sambrook J., et al. : Molecular Cloning (2nd edition) : Cold Spring Harbor Laboratory (1989).
- 2) himac APPLICATION No. 40 (1993).

HITACHI

(Export office)

Nissei Sangyo Co., Ltd.

Head Office :

1 - 24 - 14, Nishishinbashi, Minato-ku, Tokyo, 105 Japan

Tel : 81 - 3 - 3504 - 7281

Fax : 81 - 3 - 3504 - 7302

(Manufacturer)

Hitachi Koki Co., Ltd.

Scientific Instruments Div.

1060, Takeda, Hitachinaka - city, Ibaraki - pref., 312 Japan

Tel : 81 - 29 - 276 - 7384

Fax : 81 - 29 - 276 - 7475