

シールチューブ用スペースキャップ Space Caps for Seal Tubes

使用するチューブとロータの種類によりスペースキャップを選んでください。

これらの製品は、ヨーロッパ・アメリカでは、販売しておりません。

Select the correct space caps to match the tubes and rotors to be used.

These tubes and parts are not on sale in Europe and America.

チューブ Applicable tubes	スペースキャップ Space caps		適用ロータ () 内は旧形ロータ Applicable rotors (Old models)
	パーツ No. Part No.	品名 (色) Name (color)	
2PA Seal tube	S402631	Space cap (A) (blue)	(RP70VT)
5PA Seal tube	484284	B Space cap (white)	(RPV65T)
	488101	B2 Space cap (black)	P100VT, P65VT2, P90NT, P65NT2, (P83VT), (P65VF), (P55VF2), (RP65AF), (SRP83VT), (RP65VF), (RP55VF2), (RP65VT2), (RP67VF), (RP85NT)
	S403331	B3 Space cap (white)	(RP50AT4), (SRP50AT)
6.5PA Seal tube	S407828	B5 Space cap (red)	P100AT2, P50AT4, (P100AT), (RP50AT4), (RP55), (SRP50AT), (RP40-3), (RP40-2)
12PA Seal tube	484285	C Space cap (white)	P65A, (P65AT), (RP65T), (RP65), (RP40), (RPW65T), (RPW45), (P55AT)
	486582	C2 Space cap (red)	P90AT, P80AT, P70AT2, (RP83T), (RP80T-3), (SRP70AT)
	486583	C3 Space cap (white)	(RPV45T)
	S402310	C4 Space cap (black)	P65VT3, P65NT, (RP55VF), (RP65VT3), (RP65NT)
40PA Seal tube	485649	E Space cap (red)	P70AT, P50AT2, (P50A2), (P30A2), (RP70T), (RP50T-2), (RP50T), (RP50-2), (RP30-2), (RPW50T), (RP60T)
	484252	E2 Space cap (blue)	(RPV50T), (RPV30)
	S402525	E3 Space cap (white)	P50VT2, (RP50VF)
94PA Seal tube	S408969	F-Ti Space cap (red)	P45AT, (P42A), (RP45T), (RP42), (RPW35)

シールチューブを使用する際には、STF2、STF-1 形チューブシーラとチューブラックが必要です。

The tube sealer, model STF2 or model STF-1, and the tube rack are necessary to seal the tube.

Applicable seal tubes	Tube racks	Part No.	Remarks
1.5PA seal tube	Tube rack (G2)	S201778G	For STF2/STF-1
2PA seal tube	Tube rack (G)	S201778F	
3.5PA seal tube	Tube rack (B2)	S201778E	
4PA seal tube	Tube rack (B3)	S201778H	
5PA seal tube	Tube rack (B)	S201778A	
6.5PA seal tube	Tube rack (B4)	S201778J	
8PA seal tube	Tube rack (C2)	S201778L	
12PA seal tube	Tube rack (C)	S201778B	
40PA seal tube	Tube rack (E)	S201778C	
	Tube rack (F)	S201778D	
94PA seal tube	Tube rack (F2)	S201778M	For STF2
2PA,4PA,5PA,6.5PA,12PA and 40PA seal tubes	Multi-rack	S201778K	For STF2/STF-1 マルチラックでは、1.5PA、3.5PA、8PA、94PA シールチューブがご使用になれません。 The 1.5PA, 3.5PA, 8PA and 94PA seal tubes cannot be used in the multi-rack.

チューブラック G ならばラックの上面に“G”と刻印してあります。

Type of tube rack is stamped on top of the tube rack.

e.g.) In case of Tube rack (G), “G” is stamped on top of the tube rack.

